

LOUPECam®

HD CAMERAS FOR LOUPES & HEADLIGHTS

Market Leader for 10 Years • Made in USA
Internationally recognized as GoPro® Surgical • GoPro® Dental



GLASS

- HD Video 1280 x 720p
- Built-in Microphone
- Only 0.4oz/13g

Visualize • Educate • Interact • Communicate
Stream • Document • Share • Personalize • Innovate



Designs for Vision
LoupeCam®



Carl Zeiss LoupeCam®



Orasoptic LoupeCam®



Xenosys LoupeCam®



Enova XLT-125 LoupeCam®



Luxtec/Integra LoupeCam®

Technical Specifications

Weight	.46 oz / 13g
Sensor Resolution	3.0 MP
Color Depth	24-Bit True Color
Image Aspect Ratio	16:9
USB technology	2.0 High Speed
Audio	Built-in Microphone
Still Images	3.0 Mega Pixels 2048 x 1536 Pixels
Motion Video	HD 1280 x 720 Pixels@30fps
Windows File Size	2-3 GB/Hr
Mac File Size	5-6 GB/Hr



Optical Lenses

Magnification	Field of View at 16\"/>
1.5 x	10.5\"/>
2 x	5.2\"/>
3 x	3\"/>
4 x	2.5\"/>

Distributed by

System Requirements



Intel Core i5 equivalent or better
4GB RAM or more
Screen resolution of 1024 x 768 or higher
1 USB 2.0 / USB 3.0 port
1.5 GB free Hard Drive

Accessories



Bluetooth FOOT PEDAL
Hands-free camera control
and proactive editing



OPTICAL LENSES
Interchangeable
lenses for different
magnification needs



LED LIGHT
Adjustable intensity and
increased depth of field,
perfectly integrated to
work with LoupeCam®



CAMERA TRIPOD
Adaptable support
for overhead light

Mounting Adapters

Loupes

Designs for Vision
Orascopic
SurgiTel
Q-Optics
Perioptix
SheerVision
Carl Zeiss
Heine

Headlights*

Luxtec
Welch Allyn
BFW
SunOptics
Enova XLT - 125
Xenosys

(*some models, custom made
adapters also available)