

The New Generation Ligation Clip

CLIP:

Pure Titanium Clip-MRI safe

Available in 3 sizes with attractive color coding-Yellow (Extra Small), Red (small) and Blue(Medium)

DESIGN OF CLIP:

Flat rectangular body with central groove on outside which interlocks with groove on the inner side of the applicator and transverse crisscross groove on the inside -this ensures strong closure of the clip and provides a firm grip on the vessel-no slippage.

The closure is Tip to Tip.

CLIP APPLICATOR:

Color coded with thin but strong tip. For easy and smooth application especially in small and narrow access/spaces.

The inner side of clip applicator jaw is ridged which interlocks with the ridge on the outer side of the clip-thus preventing clip fall out and clip wastage.

The applicator ensures complete closure of the clip with minimal pressure applied by the surgeon.

CARTRIDGES:

Color coded with applicator.

9 clips per cartridge.

Smooth feel and easy grip.

INSTRUCTION FOR LOADING:

The applicator has to be held at 90° to the cartridge while loading the clip to avoid scissoring of the clips.

The applicator has to be held on the box like a pencil to avoid pressure on the clip while loading leading to clip loss.



Marketed by:



408-409, Orange Mall, Nr. Sharda Petrol Pump, Chandkheda,
Ahmedabad - 382424. Gujarat (INDIA)
Email: info@visionarymedtech.com • Website: www.visionarymedtech.com